

Accessories



Portable Peristaltic Pump
411720



Stainless Locking Device
310034



Wrench
310022



Flow Control Clamp
112010



Needles
14g - 113314
16g - 113316
18g - 113318
20g - 113320
22g - 113322
23g - 113323

Accessories



Sampling Kit
114018

- 1 Aseptic Septum
- 7 Syringes (60cc)
- 7 Needles (18g)
- 7 Alcohol Wipes
- 1 Dust Cover



Stainless Replacement Nut
219054



Dust Covers
310162



Quality Management, Inc.

12-Port Aseptic Fluid Access System



63-03

www.qmisystems.com

QMI Stainless Steel Fittings



Tri-Clamp End Caps (316L)

2" - 211164
2 1/2" - 211168
3" - 211174



Non-Insulated Injection Channel (316L)

212214



Insulated Tank Fitting (316L)

2" Depth - 212224
3" Depth - 212228

Custom Fittings Available Upon Request

QMI Septa

The QMI Septa are equipped with 12 individual needle guide channels. Pierce each channel once with a needle to withdraw a sterile sample. Once all channels have been used, the septa should be replaced. These septa fit inside the QMI stainless steel 12-Port fittings. Standard Operating Procedures can be found on our website at www.qmisystems.com



12 Port Septum

110064



12 Port Septum

High Temp Applications (up to 280°F)
(silicone material)

110074

QMI Collection Units



Sample Collection Bag

250ml - 111450

2L - 111730

7L - 111930

2L Peristaltic Bag 2.0mmID - 111750

2L Peristaltic Bag 3.2mmID - 111770



Syringes

60cc - 113120

140cc - 113150