

PN 110041



QUALITRU
SAMPLING SYSTEMS

INTRODUCING THE **TRU**STREAM7 SEPTUM - EU

The QualiTru TruStream7 septum EU is designed specifically for the European market, is EC1935 compliant, and is rated up to 87°C. The sterile septum has seven channels which allow for single samples to be drawn from each channel, dramatically reducing the risk of introducing contamination into the sample. Each septum features a label that clearly indicates which channels have been used. Once all channels have been used, simply replace with a new septum to continue collecting aseptic samples.



BENEFITS OF ASEPTIC INLINE PROCESS MONITORING

- Identify and mitigate in-process microbial risks early with aseptic sample points
- Monitor system hygiene through a proactive sampling and testing protocol
- Draw microbiological samples with ready-to-use sterilized equipment



**Your Test Result is Only as
Accurate as Your Sample!**

Leaders in the science of aseptic
and representative sampling.

**Call Today to Learn More or
Visit our website at QualiTru.com!**

MKTG-114 REVA 022023

Sales@QualiTru.com | 651-501-2337 | QualiTru.com