Get truly representative liquid samples as easily as 1,2,3

Sanitary Ports

QualiTru stainless steel sanitary ports are manufactured under rigorous quality control standards and allow easy access for aseptic process monitoring.





Sterile Septa

Our sterile septa allow aseptic samples to be drawn a single time from each channel virtually eliminating the risk of contamination.



We offer a complete line of single-use, sterile disposable collection units in a variety of sizes and configurations for every application.









ASEPTIC SAMPLING SYSTEM
Complete Product Guide



TRUSTREAM

TruStream7 Sanitary Ports	
TruStream7 Tank Port TruStream7 Tri-Clamp Port 1.5"	212113
TruStream7 Tri-Clamp Port 1.5"	211125
TruStream7 Tri-Clamp Port 2"	211129
TruStream7 Tri-Clamp Port 2.5"	211133
TruStream7 Tri-Clamp Port 3"	211137
TruStream7 Tri-Clamp Port 4"	211145
TruStream7 Tri-Clamp Elbow 1"	213019
TruStream7 Tri-Clamp Elbow 1.5"	
TruStream7 Tri-Clamp Elbow 2"	
TruStream7 Tri-Clamp Elbow 2.5"	213033
TruStream7 Tri-Clamp Elbow 3"	213037
TruStream7 Tri-Clamp Elbow 4"	213045
TruStream7 Recessed Port 4" x 2"	212123
TruStream7 Recessed Port 4" x 4"	212127
TruStream7 Recessed Port 5.5" x 4"	212133
TruStream7 Recessed Port 5.5" x 5"	212135
TruStream7 Recessed 5.5" x 5" w/Security Lid	212433
TruStream7 Butt Weld Elbow 1.5"	
TruStream7 Butt Weld Elbow 2"	
TruStream7 Butt Weld Elbow 2.5"	
TruStream7 Butt Weld Elbow 3"	
TruStream7 Butt Weld Elbow 4"	213095
TruStream7 Tri-Clamp Tee 1"	215141
TruStream7 Tri-Clamp Tee 1.5"	215145
TruStream7 Tri-Clamp Tee 2"	
TruStream7 Tri-Clamp Tee 3"	215151
TruStream7 Tri-Clamp Short Tee 1.5"	215515
TruStream7 Tri-Clamp Short Tee 2"	215520
To Charles 7 Double World T. 4"	245074
TruStream7 Butt Weld Tee 1"	2150/1
TruStream7 Butt Weld Tee 1.5"	
TruStream7 Butt Weld Tee 2"	
TruStream7 Butt Weld Tee 2.5"	
TruStream7 Butt Weld Tee 3"	
TruStream7 Butt Weld Tee 4"	215089
TruStream7 Wrench	310011
TruStream7 Stainless Steel Nut	
TruStream7 Dust Cover	
TruStream12 Sanitary Ports	
Trustream 17 Sanitary Ports	

2122
2111
2111
2111
2111
2122
2122
2131
2131
2131
2131
2131
2131
3100
2190
3101



TruStream Sanitary Ports

Our 3-A Certified sanitary ports are designed to work in a wide variety of locations including processing lines, silos, tanks and milk tanker trucks. QualiTru sanitary ports can also be customized to your specification.



TRUSTREAM12

QUALITRU SAMPLING SYSTEMS

Place your order today

651-501-2337

sales@qualitru.com

qualitru.com

.110011

.110021

....110084

.110064

.110074

TruStream7 Sterile Septa

TruStream12 Sterile Septa

TruStream12 Septa (high temp)....

New TruStream12 Sterile Septa.

TruStream7 Septa (high temp)...

TruStream7 Septa..

TruStream12 Septa..

TruStream 12 Channel Sterile Septa

Our 12-channel sterile septa are built for durability and easeof- use. Each single-piece septa is made of food-grade materials and features a permanent, thermal-aluminum label that clearly indicates which channels have been used. Once all 12 channels have been used, simply replace with a new sterile septa. There is no better way to gather aseptic samples.




TruStream / Channel Sterile Septa

Simple to install. Simple to use. Simple to see when a channel has been used. The TruStream7 channel sterile septa is made from food-grade contact materials and allows aseptic samples to be drawn a single time from each channel. Once all 7 channels have been used, simply replace with a new sterile septa.

QualiTru Sterile Collection Units

QualiTru offers a complete line of sterile, disposable sample bags in a variety of sizes making them ideal for all your sampling needs. All our collection units are built tough and provide a worry free handling, storing and shipping solution.



The TruDraw Single Sampler is a sterile, sealed container with an attached needle. Our new TruDraw Single Sampler is leak-proof and includes security features and ease-of-use for drawing samples up to two ounces. This industry first sampling method provides sample tracking capabilities and virtually eliminates the risk of external contamination. For added security, each TruDraw has a tamper evident feature. These simple to use samplers work with all QualiTru TruStream Septa.

TruStream sterile collection bags come in a variety of sizes and configurations for any sampling need. Each bag is equipped with a needle and tubing and



Flow Clamp Valve...

QualiTru

QualiTru Sterile Collection Units

TruStream Sampling Assembly 250ml111450

TruStream Sampling Assembly 2L111730

TruStream Sampling Assembly 7L111930

TruMotion Sampling Assembly 2L w/2.0mm silicone tubing

Single Sampler 2oz112021

for Peristaltic Pump use......111750

..112035

TruStream7 Syringe Kit #114108

TRUSTREAM7

QUALI





TruStream Collection Bags

work with all QualiTru TruStream Septa.



