

TRUST. CONFIDENCE. SAFETY.

ARE YOU QUALITRU SAFE?

THE RIGHT SYSTEMS

QualiTru Sampling Systems is committed to providing accurate, easy to use inline sampling systems for liquid Industries.

1
Port



2 Septum



3
Collection Unit



4 Pump



THE RIGHT SOLUTION,
STANDARDS AND
EQUIPMENT FOR
ACCURATE SAMPLING



QualiTru Sampling Systems offers a wide variety of products that allow safe and secure access to your fluid processing lines. These access points can be used to collect samples for a variety of applications.



TESTIMONIAL

We rely on QualiTru sampling ports throughout our process to gain better control of product quality on the front end.

”

Markets We Serve

- Craft Beer
- Dairy
- Juice
- Liquid Egg
- Milk Alternatives
- Biotech
- Nutraceuticals
- Soda
- Water
- Wine

For more information:
Sales@QualiTru.com
651-501-2337

QualiTru.com



FINALLY, GET A TRUE REPRESENTATIVE SAMPLE

That Business Decisions Can Be Made From.

OUR SYSTEM IS AS EASY AS 1-2-3-4!

1

Port



2

Septum



3

Collection Unit



4

Pump



TruStream Stainless Steel Sanitary Ports

Our ports are manufactured under rigorous quality control standards and are available in both 7- and 12-channel configurations. Custom ports can also be created for seamless integration with virtually any processing line.



TruStream Sterile Septa

QualiTru sterile septa are available in both 7- and 12-needle guide channels. Our sterile septa allow for aseptic samples to be drawn a single time from each channel. This method helps remove the risk of cross contamination due to multiple samplings from a single channel. Once all channels have been used, simply replace with a new septum.

TRUSTREAM



TruStream Sterile Collection Units

We provide a complete line of single-use, sterile, disposable collection units ranging from our 2oz TruDraw Single Sampler to our 7-Liter collection bags. Each TruStream collection bag includes tubing and needle as one unit.

TRUDRAW

TruDraw Single Sampler

The innovative TruDraw Single Sampler, with tamper-evident feature, is a sterilized, sealed container with an attached sterile needle to draw a sample from. This self-contained option virtually eliminates the risk of external contamination. The TruDraw is also ideal for simple, secure sample tracking of each individual container. Our TruDraw Single Sampler is designed to work with the TruStream sanitary ports and sterile septa.



Peristaltic Pumps

The TruMotion sampling system utilizes a peristaltic pump to draw liquid samples into a sterile TruStream collection unit.

This system is a continuous flow sampling system that collects a universal sample over the entire process run.

TruMotion collection bags for peristaltic pump use:

- 2 Liter/3.2 millimeter silicone tubing
- 2 Liter/2.0 millimeter silicone tubing

